

Quantifying uncertainties in phytoplankton  
absorption coefficients for accurate validation  
of the PACE ocean color sensor:  
moving towards satellite retrieved  
phytoplankton functional types (PFTs)





Collin Roesler

Department of Earth and Oceanographic Science

Bowdoin College

One of the PACE objectives is to support the estimation of IOPs in general and specifically optically-based Phytoplankton Functional Types




- Invert ocean color signal to retrieve hyperspectral phytoplankton absorption coefficients
- Requires capability to **validate satellite-derived phytoplankton absorption with in situ observations including uncertainties**

# Steps to construct a validation data set for $a_{\text{sat}}(\lambda)$

## 1a. Exponentially weighted $a_{\text{insitu}}(\lambda)$

- Satellite sees (spectrally) exponentially-weighted absorption



# Stepwise approach

## Step 1b. In situ hyperspectral absorption

- In situ sensors don't have the full capabilities to provide validation for PACE absorption products
  - acs - hyperspectral, but no UV and has unquantified scattering errors
  - ICAM - 9 wavelengths, 1 UV, "no" scattering error
  - a-sphere - hyperspectral but poorly constrained drift/baseline
- Validation requires merging the best of each sensor
- One approach is to quantify scattering error in acs/ac9
  - Use Integrating sphere in spectrophotometer to determine scattering correction and extrapolate to UV
  - Use ICAM to determine scattering correction and extrapolate to UV

# Objective 1b

Can the ICAM provide a pathway to expand absorption to the UV and provide scattering corrections for acs?

- Progress to date

- **Laboratory experiments** on acs and ICAM





- CDOM dilution series, additions of 50  $\mu\text{m}$  beads (Duke Scientific)
    - Phytoplankton culture dilution series

- **Field experiments** with acs and ICAM

- January 2016 Boussole
    - January 9-13 2017 closure experiments
      - Multiple acs, ac9, ICAM
      - Gershun radiometry
      - 1 psicam
      - Filter pads (Integrating sphere – 3 Pis)

# Comparison ICAM and acs







## CDOM Experiment



# Comparison ICAM and acs

## CDOM Experiment

- acs and ICAM comparison





Wavelength (nm)

# Comparison ICAM and a<sub>cs</sub>

## CDOM Experiment



- Water calibrated, a<sub>cs</sub> temperature corrected
  - Comparable values
  - Scales with bead additions





# Comparison ICAM and acs

## CDOM Experiment

- Water calibrated, acs temperature corrected
  - Comparable values
  - Scales with bead additions






- acs scattering corrected
  - $a_{ICAM} > a_{acs}$
  - acs independent of bead addition





# Comparison ICAM and acs

## Phytoplankton Experiment



# Comparison ICAM and acs Phytoplankton Experiment

- acs and ICAM comparison with phytoplankton culture
  - ICAM overestimates acs by 5 to 14%  $\pm$  3 to 11% ( $r^2 > 0.99$ ,  $p < 10^{-6}$ ) at all wavelengths
  - But underestimates acs(676) by 55%  $\pm$  3% ( $r^2 > 0.99$ ,  $p < 10^{-6}$ )
  - Is fluorescence causing the underestimation?



# Summary of acs and ICAM Comparison

- Wavelength range and resolution
  - acs ~80 channels, 400 to 750 nm (3 nm but requires unsmoothing)
  - ICAM 9 channels, 365 to 676 nm (bandwidth not reported)
- Data processing and corrections
  - Pure water correction
  - Temperature and salinity correction (acs published protocol)
  - Scattering correction (acs uses beam attenuation)
- Data rates
  - acs 6 Hz
  - ICAM 1 Hz (now 0.5 to 0.2 Hz)
- Uncertainties
  - ICAM exhibits sensitivity to particle scattering but does not have sufficient wavelength range for developing even a null correction minimally
  - ICAM 676 nm channel significantly underestimates the red chlorophyll absorption peak by >50%
- In situ closure experiment conducted last week will investigate whether these are systematic issues

# Stepwise approach

## Step 2. Discrete analysis of phytoplankton absorption

- Samples are concentrated on glass fiber filters
  - Increases signal to noise
  - Removes CDOM from analysis
  - Separation of phytoplankton from other particulates
  - But there are scattering losses and pathlength amplification
- Improved spectrophotometric configuration
  - Samples center-mounted in integrating sphere (IS)
  - Mitigates scattering losses
  - Reduces uncertainty in pathlength amplification

**Objective: Quantify uncertainties in estimating phytoplankton absorption coefficients on discrete samples**

# Objective 2





Quantify uncertainties in estimating phytoplankton absorption coefficients on discrete samples

## Progress to date

- Transfer function between transmission mode and integrating sphere measurements




- Reduce uncertainty and increase accuracy of historical filter pad absorption data

- Range of filter loading
- Field/lab samples



# Transfer function between transmission mode and integrating sphere measurements


- (left) Results for 30 samples ( $n > 14,000$ ) over the entire spectrum
- (right) restricting observations  $OD < 0.35$  (%T mode)
- dotted is the best fit linear model  
 $OD_{IS} = 1.29 * OD_{TR} + 0.00205$
- Tested with data from SOCLIM (La Réunion-Kerguelen 10/2016)




# Objective 3

## Discerning absorption-based phytoplankton functional types from multispectral fluorescence

- Some Validation Platforms (floats, moorings) not optimal for ac meter deployments (e.g. size, power, data, biofouling)
- To increase opportunities for validation, small, low power optical sensors are **also** needed
- Multichannel fluorometers (3X1M) provide means for discerning absorption ratios and pigmentation (*Proctor and Roesler 2010; Thibodeau et al. 2014*)




← Fluorescence is proxy for absorption  
→ Fluorescence ratios are proxies for pigment composition



# Objective 3

## Discerning absorption-based phytoplankton functional types from multispectral fluorescence


- Problems linking in situ fluorescence observations with remote sensing observations due to Non-Photochemical Quenching (NPQ)
- Remote sensing requires daytime validation and unquenched fluorescence occurs at nighttime
- Particularly problematic for tidally impacted coastal environments



# Objective 3

Discerning absorption-based phytoplankton functional types from multispectral fluorescence


- Start with hourly fluorescence observations



# Objective 3

Discerning absorption-based phytoplankton functional types from multispectral fluorescence


- Remove NPQ-impacted observations (grey)



# Objective 3

Discerning absorption-based phytoplankton functional types from multispectral fluorescence


- Use tide proxy (ADCP currents or Salinity) to identify hours of high (blue) and low (red) tide *during non-NPQ times*



# Objective 3

Discerning absorption-based phytoplankton functional types from multispectral fluorescence


- Interpolate daily tidal excursions in non-NPQ fluorescence between high and low tide envelope to estimate the non-quenched calibrated fluorescence at the mooring that would be observed every hour



# Objective 3

Discerning absorption-based phytoplankton functional types from multispectral fluorescence


- Improved matchup between hourly in situ calibrated chlorophyll fluorescence time series (black line) and radiometric chlorophyll concentration (green symbols)



# Objective 3

Discerning absorption-based phytoplankton functional types from multispectral fluorescence


- Now variability between  $\text{Chl}_{\text{sat}}$  and  $\text{Chl}_{\text{insitu}}$  can be examined from the perspective of ***natural variability in fluorescence per chl***



# Objective 3

## Discerning absorption-based phytoplankton functional types from multispectral fluorescence

- Natural variations in the fluorescence per chl (calibration slope) are well documented
- WETLabs calibration overestimates in situ chl with a global bias of 2 but local biases of over 6 in some regions



Challenges of and recommendations for absolute calibration of in situ fluorometers for assessing global phytoplankton distributions

Roesler et al 2017 accepted to Limnology and Oceanography Methods

Thank you

